

CASE STUDY

How to surpass significantly the high demands of ISO 8968-1 for protein with a simple Kjeldahl method



Abstract

With the target to optimise the range of variation in Kjeldahl protein analysis of milk according to official norms like AOAC 991.20 or ISO 8968-1:2001, digestion conditions were varied to obtain lowest variation in protein content between duplicates (significantly below 0,038 % protein) with highest possible recovery .

It has been found that digesting for 2 hours at 390 °C with a KJELDATHERM® digestion block KBL20s gives very convincing results, especially in the range of variation the maximum allowed value of 0,038 % protein between duplicates was significantly lower, the maximum range could be even kept in series of 20(!) samples.

The only difference found compared to the mentioned official norms is the status of the sample after digestion. Although the norms describe a clear solution, we could obtain very good results with the parameters mentioned below, even with samples which were still dark after digestion. The quality of the results justify this small difference.

Recovery was checked with L-Tryptophane (AppliChem, purity min. 99 %). Average recovery found was 99,494 % with the described method below which is significantly higher as 98 % as demanded for L-Tryptophane by the norms.

Method Data and Results

Sample type:	Milk
Submission no.:	3590, 3667
Date:	20.09.-19.10.2012

Analytical details:

Chemicals

Water:	demineralised or distilled
Sulphuric acid:	H ₂ SO ₄ 98 % at least
Catalyst tablets:	Kjeltabs CX analog Kjelcat Cu
Caustic-Soda:	NaOH 32 %
Boric acid:	H ₃ BO ₃ 2 %

Standard Solution:

Hydrochloric acid:	HCl c=0.1 mol/l or
Sulphuric acid:	H ₂ SO ₄ c=0.05 mol/l

All chemicals applied have p.a. grade.

Instruments

Analytical Balance	(0.0001 g)
Kjeldahl Digestion Block	KJELDATHERM®
KBL20s, TURBOSOG-Scrubber	
Distillation system	VAPODEST® 50 s c

Sample preparation

Heat up the milk in a water bath of 40 °C to a temperature of 38 °C, stir up carefully and let cool down to room temperature. The initial sample weight of liquid samples is determined with a disposable syringe by differential weighing. Sulphuric acid is used to wash down any sample residue, which might remain at the glass walls.

Sample:	~4g
Catalyst:	2 Kjeltabs type CX or analog 2 Kjelcat Cu
H ₂ SO ₄ :	20 ml min. 98 %



Method	Series Position	Date	Sample	Weight [g]	Blank value [mL]	Consumption [mL]	N [%]	P [%]	Duplicate determination	
Milk	1	18.10.2012 12:11	3667	3,4133	0,2023	12,8085	0,5173	3,3005		
Milk	2	18.10.2012 12:25	3667	3,6963	0,2023	13,8524	0,5173	3,3002	0,0003	
Milk	3	18.10.2012 12:32	3667	3,5943	0,2023	13,5570	0,5204	3,3204		-0,0202
Milk	4	18.10.2012 12:39	3667	3,8219	0,2023	14,2929	0,5164	3,2947	0,0257	
Milk	5	18.10.2012 12:46	3667	3,6595	0,2023	13,7306	0,5178	3,3036		-0,0089
Milk	6	18.10.2012 12:53	3667	3,7314	0,2023	14,0334	0,5192	3,3125	-0,0089	
Milk	7	18.10.2012 13:00	3667	3,5626	0,2023	13,3612	0,5174	3,3008		0,0117
Milk	8	18.10.2012 13:07	3667	3,8268	0,2023	14,3613	0,5183	3,3065	-0,0056	
Milk	9	18.10.2012 13:15	3667	3,8050	0,2023	14,2601	0,5175	3,3016		0,0048
Milk	10	18.10.2012 13:21	3667	4,4999	0,2023	16,8589	0,5185	3,3079	-0,0062	
Milk	11	18.10.2012 13:28	3667	4,0541	0,2023	15,2076	0,5184	3,3076		0,0003
Milk	12	18.10.2012 13:35	3667	3,9354	0,2023	14,8308	0,5207	3,3218	-0,0142	
Milk	13	18.10.2012 13:42	3667	3,5699	0,2023	13,4644	0,5204	3,3199		0,0019
Milk	14	18.10.2012 13:48	3667	3,9336	0,2023	14,7559	0,5182	3,3063	0,0135	
Milk	15	18.10.2012 13:55	3667	3,9511	0,2023	14,8226	0,5183	3,3068		-0,0004
Milk	16	18.10.2012 14:02	3667	3,8191	0,2023	14,3456	0,5187	3,3095	-0,0027	
Milk	17	18.10.2012 14:09	3667	3,9122	0,2023	14,6109	0,5159	3,2913		0,0182
Milk	18	18.10.2012 14:16	3667	3,8360	0,2023	14,4045	0,5186	3,3086	-0,0173	
Milk	19	18.10.2012 14:22	3667	3,8617	0,2023	14,5230	0,5194	3,3140		-0,0054
Milk	20	18.10.2012 14:29	3667	3,9812	0,2023	14,9423	0,5186	3,3086	0,0054	
						Maximum value	0,5207	3,3218		
						Average value	0,5184	3,3071		
						Minimum value	0,5159	3,2913		
						Range of Variation	0,0048	0,0305		

Method	Series Position	Date	Sample	Weight [g]	Blank value [mL]	Consumption [mL]	N [%]	P [%]	Duplicate determination	
Milk	1	19.10.2012 12:13	3667	3,7384	0,2023	14,0532	0,5190	3,3110		
Milk	2	19.10.2012 12:19	3667	3,9010	0,2023	14,7057	0,5208	3,3225	-0,0115	
Milk	3	19.10.2012 12:26	3667	3,8669	0,2023	14,5834	0,5209	3,3235		-0,0010
Milk	4	19.10.2012 12:33	3667	3,8625	0,2023	14,5464	0,5202	3,3187	0,0048	
Milk	5	19.10.2012 12:40	3667	3,9185	0,2023	14,7356	0,5195	3,3144		0,0043
Milk	6	19.10.2012 12:47	3667	3,7295	0,2023	14,0657	0,5207	3,3219	-0,0075	
Milk	7	19.10.2012 12:53	3667	4,1781	0,2023	15,7384	0,5208	3,3230		-0,0011
Milk	8	19.10.2012 13:00	3667	4,0103	0,2023	15,1072	0,5206	3,3214	0,0016	
Milk	9	19.10.2012 13:07	3667	4,0094	0,2023	15,0517	0,5188	3,3097		0,0116
Milk	10	19.10.2012 13:14	3667	3,9677	0,2023	14,8945	0,5187	3,3091	0,0006	
Milk	11	19.10.2012 13:21	3667	3,9054	0,2023	14,7275	0,5210	3,3237		-0,0146
Milk	12	19.10.2012 13:28	3667	4,1788	0,2023	15,7528	0,5212	3,3255	-0,0018	
Milk	13	19.10.2012 13:34	3667	4,1184	0,2023	15,5028	0,5204	3,3200		0,0055
Milk	14	19.10.2012 13:59	3667	4,0040	0,2023	15,0695	0,5201	3,3182	0,0018	
Milk	15	19.10.2012 14:06	3667	4,0282	0,2023	15,1327	0,5192	3,3123		0,0059
Milk	16	19.10.2012 14:13	3667	4,1692	0,2023	15,6593	0,5193	3,3131	-0,0009	
Milk	17	19.10.2012 14:20	3667	3,9925	0,2023	15,0162	0,5197	3,3158		-0,0027
Milk	18	19.10.2012 14:27	3667	4,1946	0,2023	15,8242	0,5217	3,3282	-0,0124	
Milk	19	19.10.2012 14:33	3667	3,9430	0,2023	14,9214	0,5229	3,3360		-0,0078
Milk	20	19.10.2012 14:40	3667	4,2586	0,2023	16,0484	0,5212	3,3252	0,0107	
						Maximum value	0,5229	3,3360		
						Average value	0,5203	3,3197		
						Minimum value	0,5187	3,3091		
						Range of Variation	0,0042	0,0268		



C. Gerhardt – Quality made in Germany

AUTOMATING STANDARD ANALYSES

Completely automated laboratory analysis systems from C. Gerhardt are highly developed special equipment. They automate recurring analysis processes in accordance with national and international standards and norms. They continuously provide precise and reproducible analysis results quickly, at low cost, economically and highly efficiently.



An excerpt from our product portfolio:

- + **COMPLETELY AUTOMATIC HYDROLYSIS**
HYDROTHERM – automatic acid hydrolysis system for fat determination according to Weibull-Stoldt. When combined with SOXTHERM, HYDROTHERM is an ideal system solution for total fat determination.
- + **COMPLETELY AUTOMATIC FAT EXTRACTION**
SOXTHERM® – automatic fast extraction system for fat determination
- + **COMPLETELY AUTOMATIC CRUDE FIBRE EXTRACTION**
FIBRETHERM® – completely automated processing of the boiling and filtration processes for determining crude fibre, ADF and NDF.
- + **COMPLETELY AUTOMATIC NITROGEN ANALYSIS**
DUMATHERM® – nitrogen/protein determination of solid and liquid samples according to the Dumas combustion method. A fast and convenient alternative to the classic Kjeldahl method for almost all sample matrices.

C. Gerhardt GmbH & Co. KG

Cäsariusstraße 97 • 53639 Königswinter, Germany
Tel.: +49 (0) 2223 2999-0 • www.gerhardt.de

Presented by:

