

Exhaust equipment

TURBOSOG



Application

Centrifugal scrubber for the condensation and neutralization of aggressive acid fumes. The TURBOSOG works in two steps, separating and washing out acid fumes. The removal of the fumes is very efficient with extremely low running costs. Little service is required as no activated carbon filters are used in the system.

- TURBOSOG scrubber unit, order no.12-005

Product features

- compact bench top unit
- extremely low running costs
- very efficient fume removal
- virtually maintenance-free

Recommended accessory for digestion of aqueous samples

Additional condensate bottle for digestions of aqueous samples with **KJELDATHERM** or **TURBOTHERM** systems. Easy handling with screw caps and snap lids. The unit is attached to the side of the TURBOSOG and connected to the water supply. Existing units can be upgraded at any time.

- Additional cooling unit for TURBOSOG, order no. 12-0361

