

Traditional Kjeldahl distillation apparatus



Kjeldahl distillation apparatus available with 50 - 750 ml Kjeldahl flasks.

Distillation apparatus with Kjeldahl flasks

50-250 ml

Serial heating unit KI 16 for 6 distillations in Kjeldahl flasks 50 ml or 100 ml (micro) resp. 250 ml (semi micro). Supplied with Kjeldahl flasks, Reitmair tops, condensers, outlet tubes, Erlenmeyer flasks, all rubber connections, support rings and power cord.

Available models

- KI 9/16 Distillation apparatus with 6 flasks 50 ml, order no. 10-0037
- KI 9/16 Distillation apparatus with 6 flasks 100 ml, order no. 10-0038
- KI 9/16 Distillation apparatus with 6 flasks 250 ml, order no. 10-0039

Distillation apparatus with Kjeldahl flasks 500-750 ml

Serial heating unit KI 24 or KI 26 for 4 resp. 6 distillations in Kjeldahl flasks 500 ml or 750 ml (macro). Supplied with Kjeldahl flasks, Reitmair tops, condensers, outlet tubes, Erlenmeyer flasks, all rubber connections, support rings and power cord.

Available models

- KI12/24 Distillation apparatus with 4 flasks 500 ml, order no. 10-0040
- KI12/24 Distillation apparatus with 4 flasks 750 ml, order no. 10-0041
- KI12/26 Distillation apparatus with 6 flasks 500 ml, order no. 10-0042
- KI12/26 Distillation apparatus with 6 flasks 750 ml, order no. 10-0043

[Click here](#) for information on automatic distillation units: VAPODEST.