

Software for rapid extraction systems

SOX THERM Manager



SOX THERM-Manager is a powerful, multitasking software, which meets all the demands of a modern laboratory management.

The installation on every common PC or Laptop under Microsoft Windows 98, 2000 or XP is easy and fast.

In addition to the functions of the MULTISTAT controlling unit, multiple documentation and logging options, like e.g. the sample name, program and recognition of the user are offered as well.

The functions

- Control of up to four SOX THERM units (up to 24 simultaneous extractions), each unit with individual programming
- Monitoring of the SOX THERM units
- Programming and storage of the extraction programs
- Logging and documentation of all relevant data
- Diagnosis function for an ideal extraction process
- Concise arrangement of all relevant status displays in separate windows facilitates the best possible supervision

Everything under Control

■ Clearly structured program windows

Each SOX THERM unit has a program window and a status window. At any given time, the user has a perfect overview of all functions and programs in use. Up to 4 SOX THERM units are displayed on one screen.

■ Start immediately

You are working with a windows user interface, so that the program looks familiar right away. This saves time-consuming familiarization periods.

■ Individual settings

It is up to you to determine, which functions should be displayed permanently in the status window. These windows can be displayed side by side, or one after the other, can be shown or hidden.

■ Program library

For each extraction unit you can start an individual program library. Upon request, we will send you application notes for all established procedures.

■ Documentation

Thanks to vast logging functions, the SOX THERM Manager offers a vast variety of possibilities for documentations, like e.g. starting time, operator name, entire running time, error messages, comments and many more.