

System for the determination of ADF, NDF and crude fibre

FIBRE THERM



System developed to automate the time consuming dosing, boiling and filtration steps required in crude fibre, ADF and NDF analysis.

System configuration FIBRE THERM

- Closed system
- 9 programs
- Automatic acid addition
- Automatic caustic addition
- Automatic rinse cycle
- Automatic suction
- Cooling water control
- Complete function control
- Batch handling of 12 samples

Optional accessories for FIBRE THERM

- External automatic amylase dosing pump (required in NDF determinations)
- Set of tanks with level sensor control

Imperative accessories for FIBRE THERM

- 100 pieces FIBREBAGS for crude fibre determination, order no. 10-0128
- 100 pieces FIBREBAGS for ADF and NDF determinations, order no. 10-00127
- Glass spacer, order no. 10-0124