

Traditional apparatus for crude fibre determination



For the determination of crude fibre in feed.

Complete unit ready to operate with special supports, glass condensers, ring shaped rests, silicone rubber tubings, beakers 1000 ml and power cord.

Available models

- RF 16/1 Crude fibre digester, 1 place, 230 V 50-60 Hz, 450 W, order no. 10-0016
- RF 16/4 Crude fibre digester, 4 places, 230 V 50-60 Hz, 1800 W, order no. 10-0017
- RF 16/6 Crude fibre digester, 6 place, 230 V 50-60 Hz, 2700 W, order no. 10-0018

Information on systems with programmable dosing, boiling and filtration steps can be obtained here: [FIBRE THERM](#).