

## Hotplates

EV



Single and serial heaters with hotplates of 85 mm diameter. Each hotplate with graduated fully adjustable energy controller and excess current cut off with pilot lamp. Maximum temperature 425 °C.

### Application

Thanks to the broad variety of accessories the hotplates can be used as:

- Extraction units
- Crude fiber units
- Hydrolysis units
- Digestion units
- Distillation units

### Available models

- EV1 Single hotplate, 1 place, 230 V 50-60 Hz 450 W, order no. 10-0001
- EV14 Serial hotplate, 4 places, 230 V 50-60 Hz 1800 W, order no. 10-0005
- EV16 Serial hotplate, 6 places, 230 V 50-60 Hz 2700 W, order no. 10-0006
- other voltages **on request**



### Accessories for the heating of beakers and flasks with flat bottom

- Top mould diameter 120 mm for beakers up to 2000 ml and flasks with flat bottom, order no. 10-0051

### Accessories for heating of flasks

- Top mould diameter 65 mm for flasks up to 100 ml, order no. 10-0052
- Top mould diameter 95 mm for flasks 250 ml to 500 ml and Kjeldahl flasks 500 ml to 750 ml, order no. 10-0053
- Airbath insert, order no. 10-0056
- Support rod 600 mm, order no. 10-0097
- Support rod 1000 mm, order no. 10-0100
- Forked rod with screw thread, order no. 10-0098