

RAPID DRYER FOR LABORATORY GLASSWARE



Laboratory glass dryer in tried and tested C. Gerhardt quality.

Dry laboratory glassware in just 6 – 15 minutes thanks to filtered hot air.

Simultaneous drying of up to 56 tubes.
Rotatable dryer attachment with 24 air nozzles with 16 mm diameter and 32 nozzles with 11 mm diameter.
Maximum temperature in tube approx. 50 °C.

Cold or hot operation. With electrical timer.

Scope of delivery Rapid glass dryer
Instruction manual

Order information

Order No.	Type	Description
10-0010	STL56	Rapid Glassdryer, 230V, 50Hz Version 115 V
10-0200	STL56	Rapid Glassdryer, 115 V, 60 Hz

Technical data

Type	STL56
Electric Timer [min]	0 - 120
Drying tubes	56
Dimensions W / D / H [mm]	380 / 380 / 1150
Weight [kg]	12
Nominal voltage [V]	230 / 115
Frequency [Hz]	50 / 60
Nominal wattage [W]	1000



C. Gerhardt GmbH & Co. KG
Cäsariusstraße 97 · 53639 Königswinter, Germany
Tel.: +49 (0) 2223 2999-0 · www.gerhardt.de

