


High-quality consumables



The performance of an analytical system is only as good as its consumables. That is why our consumables have been specifically developed to complement our analytical systems. Because precise analytical results can only be guaranteed if consumables and systems are designed to work together.

Fat determination and extraction

Consumables for your daily work:



Folded filters

- + Specially coated to capture hydrolysed fat quantitatively
- + Time saving through efficient sample processing thanks to ideal pore size
- + Specially designed for use in HYDROTHERM

Order number 1004092



Weighing paper

- + Practical block for tearing single papers off
- + Easy to shape after sample weight
- + Almost residue-free decomposition
- + Free of grease
- + Protects the sample from burning and ensures a slow decomposition for a controlled reaction

Order number 1004939






Extraction beakers

- + Robust quartz glass and optimum glass thickness result in a long life time and consistent boiling process
- + Fire-polished beakers ensure reliable sealing for less evaporation of solvent
- + Three different options: stainless steel, brown glass, quartz glass tailored to your application


Order number 13-0050 (makro)
Order number 13-0077 (makro with barcode)
Order number 13-0080
(brown, light-sensitive samples)
Order number 13-0051
(mikro, reduced use of solvent)



Extraction thimbles

- + Reusable
- + Flexibility thanks to different sizes
- + Optimal pore size for reliable solvent flow
- + Durable glass sleeves are acid and heat resistant
- + Ideal positioning of thimble using special holder


Order number 13-0054
(made of cellulose, 33 x 80 mm)
Order number 13-0057
(made of cellulose, 33 x 94 mm)
Order number 13-0095
(made of glass, 34 x 94 mm)
Order number 13-0096
(made of glass, 40 x 94 mm)



Holder for extraction thimble

- + Specially designed for SOXTHERM beakers
- + Optimisation of the solvent flow through ideal positioning in the extraction beaker

Order number 13-0061 (for thimbles Ø 25 mm)
Order number 13-0062 (for thimbles Ø 33 mm)
Order number 13-0063 (for thimbles Ø 40 mm)



Boiling stones

- + Large surface area prevents delayed boiling and foaming
- + No decomposition or reaction with the solvent

Order number 1000774 (Typ A, 250 g)
Order number 13-0047 (resistent, 50 g)



Consumables by CG

- Top-quality analyses
- High-quality reagents
- Guaranteed performance with C. Gerhardt systems
- In accordance with ISO 9001
- Proven development process
- In-house quality assurance



Reliable partner

- Innovative products
- Many years of experience
- Personalised consultation
- Training for your employees
- Attractive quantity discounts



Our promise to you

- Smooth order process
- Secure supply chain
- Reliable availability
- Flexible service concepts
- Made in Europe
- Excellent price-performance ratio

You can also find all of our consumables on our website. Conveniently put together a selection of consumables for your daily work here.



www.gerhardt.de/en/products/accessories-consumables/

Contact:
lead@gerhardt.de

or your authorized local distributor

www.gerhardt.de/en