

High-quality consumables



The performance of an analytical system is only as good as its consumables. That is why our consumables have been specifically developed to complement our analytical systems. Because precise analytical results can only be guaranteed if consumables and systems are designed to work together.

Reliable fibre analysis

Consumables for your daily work:



FibreBag

- + Standardisation of filtration through defined mesh size
- + Large surface area for effective detergent flow without clogging
- + Residue-free incineration
- + Sample is permanently rinsed
- + Nitrogen-free, making it suitable for determining the content of cell wall-bound proteins (e.g. NDLP), for protein determination using Kjeldahl in the context of fibre determination and for chitin analysis in insect-based animal feeds
- + Adapted pore size tailored to the demands of CF, NDF, ADF or ADL determination

Order number 10-0127 FibreBags ADF, for ADF-/NDF determination

Order number 10-0128 FibreBags RF, for crude fibre determination

Order number 10-0142 FibreBags S, for very fine grinded samples





Quartz glass crucible

- + Visual control of the analytical progress
- + Optimal size for an analysis
- + Long durability due to robust quartz glass
- + Optimisation of the workflow


Order number 13-0048



Glass spacer

- + Optimal rinsing of the sample
- + Optimised sample weight in the FibreBags
- + Long life time due to robust glass


Order number 10-0124



Anti-foam solution

- + Prevents the sample from foaming
- + Ready for immediate use
- + Simple dosing because of dosing bottle

Order number 10-0143



Reference substance

- + With certificate
- + To qualify your analytical results

Order number 10-0129 reference substance 1:
for determining protein, fat and crude fibre, 30 g

Order number 10-0130 reference substance 2:
for determining protein, fat and crude fibre,
ADF/NDF, 30g



Consumables by CG

- Top-quality analyses
- High-quality reagents
- Guaranteed performance with C. Gerhardt systems
- In accordance with ISO 9001
- Proven development process
- In-house quality assurance



Reliable partner

- Innovative products
- Many years of experience
- Personalised consultation
- Training for your employees
- Attractive quantity discounts



Our promise to you

- Smooth order process
- Secure supply chain
- Reliable availability
- Flexible service concepts
- Made in Europe
- Excellent price-performance ratio

You can also find all of our consumables on our website. Conveniently put together a selection of consumables for your daily work here.



[www.gerhardt.de/en/products/
accessories-consumables/](http://www.gerhardt.de/en/products/accessories-consumables/)

Contact:
lead@gerhardt.de

or your authorized local distributor

www.gerhardt.de/en